

Our Products

Catheters



Hematuria Catheter



Intended Use :

hematuria drainage of the bladder or post-operative drainage and irrigation of the bladder after prostatectomy.

Size

6Fr, 8Fr, 10Fr, 12Fr, 14Fr, 16Fr
18Fr, 20Fr, 22Fr, 24Fr, 26Fr

Ingredients

High Grade Silicone

Key Features

- 100% Silicone Material
- Radiopaque (X-ray) Line
- MRI-Safe Valve
- Enhanced Anti-Blockage Drainage Lumen

