

Our Products

Stents



Key Features

- Thermosensitive, advanced PU material for enhanced patient comfort
- Stent will be stiffer at room temperature promoting pushability for stent deployment
- Material becomes softer at body temperature
- Optimal loop shape and high memory retention reduces the risk of migration
- Optimal radiopacity for better visualization

Standard	
Size (in Fr)	Length (in cm)
4.5, 4.7, 4.8, 5, 5.5, 6, 6.5, 7 & 8	12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

Multiloop	
Size (in Fr)	Length (in cm)
4.7, 4.8, 5, 6, 7 & 8	18-24 & 24-30

Pediatric	
Size (in Fr)	Length (in cm)
3, 3.5, 3.7, 3.8 & 4	6, 8, 10, 12, 14 & 16