

Our Products

Stone Basket



Delta Basket



Intended Use :

Steerable Basket for rigid nephroscope for stone extraction from calices which are inaccessible due to acute angle.

Size (in inch)

7

Length (in cm)

40

Design Features

- Can be angled inside kidney to reach inner calyx
- Can be used for straight path stone manipulation and also can be angled both direction
- Unique basket shape which prevents stone migration

Procedural Features:

- Prevents injury to the infundibulum and reduce bleeding
- Avoids the need for multiple tracts
- Can be used without the assistance of flexible nephroscope
- Helps to remove big blood clots