

Our Products

Hemostats



FAM Trium

Advanced multilayer haemostat



Intended Use :

Controls capillary, venous, and small arterial bleeding; prevents early postoperative bleeding in parenchymal, muscular, and cavity surgeries, including classical, robotic, and endoscopic procedures.

Size

2.5cm X 5cm, 5cm X 7.5cm
5cm X 10cm, 10cm X 10cm

Ingredients

Oxidised Non Regenerated Cellulose

- Key Features**
- Multilayer Fibrous Structure
 - Fast and Effective Hemostasis
 - High Blood Absorption
 - Fully Resorbable Material
 - Bactericidal & Antifungal Action
 - Flexible & Trocar-Friendly Design

