

Our Products



Guidewires

PTFE Coated Coil



Types

J Tip | Straight Tip

Size (in inch)

0.018 | 0.021 | 0.025 | 0.032 | 0.035 | 0.038

Length (in cm)

150 | 175 | 260

Intended Use :

Provide a smooth, safe, and supportive pathway for the insertion and navigation of catheters or other medical devices through vessels or anatomical channels during diagnostic or therapeutic procedures.

Key Features :

- Stainless steel core construction
- New thumb indentation near Luer opening
- Atraumatic "J" tip design
- Ultra-soft tip with ultra-short flexibility
- Radiopaque tip for precise visibility.

Intended Use :

To guide and support catheters or other devices into blood vessels during procedures and to navigate difficult, narrow, or tortuous vessels with low friction

Key Features :

- Hydrophilic coating for excellent lubricity.
- 7cm floppy tip for easy, atraumatic insertion.
- Excellent torqueability for precise control.
- Nitinol core for kink resistance.

Hydrophilic Guidewire



Types

J Tip | Straight Tip | Angeled

Size (in inch)

0.032 | 0.035 | 0.038

Length (in cm)

150 | 260