


Our Products

Catheters



Hemodialysis Catheter



Intended Use :

Designed for short-term vascular access during dialysis and apheresis therapy, for intravenous or percutaneous insertion into the patient's jugular, subclavian, or femoral vein.

Size:

12Fr x 15cm, 12Fr x 16cm, 12Fr x 20cm, 12Fr x 24cm

Type:

Curved, Straight

- Key Features**
- Medical-grade Polyurethane Catheter
 - Biocompatible and Hemocompatible
 - Nickel Titanium Alloy Guidewire

