


Our Products

Interventional devices



Introducer Sheath

PEELABLE SHEATHS WITH DILATOR



Intended Use :

A peelable sheath with dilator is used to gently dilate the entry site and provide a stable access pathway for inserting medical devices, and it can be easily peeled away without disturbing the device after placement.

Key Features :

- Smooth insertion
- Reduced tissue trauma
- Peel-away design
- Medical-grade materials
- Safe and reliable

Size (in inch)

Sheath S & L: 4Fr X 7cm, 5Fr X 7cm, 6Fr X 7cm, 7Fr X 7cm

In Package

All include J Tip (0.021" X 45CM)

Introducer needle= 20G X 4CM