


Our Products

Interventional devices



Introducer Sheath

RADIAL SHEATHS WITH DILATOR



Intended Use :

To provide a stable, low-trauma entry pathway into a blood vessel, enabling safe insertion and guidance of catheters, guidewires, or interventional devices during diagnostic or therapeutic procedures.

Key Features :

- Safe Vascular Access
- Hemostasis Valve
- Compatibility & Flexibility
- Smooth Insertion
- Radiopaque Markings

Size (in inch)

Sheath S & L: 5Fr X 7cm, 6Fr X 7cm, 7Fr X 7cm

In Package

All include J Tip (0.025" X 45CM)

IV Catheter 20G X 4.5CM