

Our Products

Leukocyte Reduction Filter



Intended Use :

For the removal of leukocytes (white blood cells) from blood and blood components during transfusion or blood processing to reduce the residual leukocyte count, minimizing transfusion-related complications associated with donor leukocytes.

Type:

- Leukocyte reduction filter for RBC (bedside)
- Leukocyte reduction filter for platelet (bedside)
- Leukocyte reduction filter for RBC (blood bank)

Key Features

- Easy handling, no priming required
- High filtration performance, 100% microaggregates removal
- Better blood compatibility, gamma rays and EO sterilised

