


Our Products

Polypropylene Mesh



Intended Use :

Surgical mesh (also called hernia mesh) is a medical implant made of loosely woven mesh, which is used in surgery as either a permanent or temporary structural support for organs and other tissues

Size

6x11cm, 7.5x15cm, 10x15cm, 15x15cm, 15x20cm, 15x30cm, 20x20cm, 20x30cm, 30x30cm

Ingredients

Polypropylene

Key Features

- Biologically compliant.
- Flexible, robust and lucid..
- High power and displosion resistance for sustainable support.
- Does not cluster and stays strong for long periods of time.

