

Our Products

Soft Tissue Biopsy Needles



Intended Use :

For the percutaneous core tissue specimens from soft tissue lesions or organs for histological examination, to obtain intact tissue samples in the diagnosis of benign or malignant diseases.

Size (working length):

14G, 16G, 18G, 20G x
10cm, 15cm, 20cm

Type:

Without coaxial needle,
with coaxial needle

- Key Features**
- Choice between two different sampling depths
 - Atraumatic sharpened tip
 - Slide stopper
 - Internal echo-marker depth marks

