


Our Products

Catheters



Ureteral Catheter



Intended Use :

Used for catheterization of the urinary tract during intraoperative procedures. Used for injection of contrast media.

Size (in inch)

3, 4, 5, 6, 7 & 8

Length (in cm)

70

Key Features

- Optimal radiopacity for better visualization
- With 1 cm marking which helps in catheter position, measurement of ureter length
- Detachable hub for scope insertion

Type

- Straight Tip Open End
- Straight Tip Closed End
- Cone Tip Open End
- Cone Tip Closed End
- Angled Tip Open End
- Angled Tip Closed End